

Study Design:

Verify that the candidate reference interval is in concordance with ASVCP guidelines. Verify that the patient population, preanalytical and analytical factors from the candidate reference interval are the same or very similar to the receiving laboratory.

Step One

Collect 20 results from the laboratory's healthy patient population (n=20) or reference individuals (n=20)

Step Two

Assess for outliers of patient results using Dixon's or Tukey's outlier analysis. Discard outliers and replace results.

Step Three

Assess the distribution and count the number of patient results that fall outside the candidate reference interval

Results are distributed across the range of the candidate RI AND 1 or 2 results fall outside of the candidate RI

Distribution is too narrow or patient results are clustered toward one end of the candidate RI OR zero or ≥ 3 results fall outside of the candidate RI
Transference NOT verified

Reject candidate RI

Step Four

Collect 20 additional samples from a new cohort of reference individuals

If 1 or 2 results fall outside of the candidate RI

If zero or ≥ 3 results fall outside of the candidate RI
Transference NOT verified

These 40 patient results can be used to construct a denovo RI (n= 80 is recommended)

Verified!