



www.immunizationmanagers.org

Vendor: LogTag/ Control Solutions		
Product Type: Vaccine Temperature Monitoring Data Logger		
Product Name: VFC 400 Vaccine Temperature Monitoring Kit		
Vendor Primary Presenter: Steve Wick Control Solutions		
Corporate Alliance Member: LogTag Recorders		
Immunization Programs that currently use product: 36		
CDC Recommendations		
	Yes	No
Meets or exceeds all CDC guidelines for continuous temperature monitoring	Y	
Digital display on the outside of the storage unit to allow users to read temperatures without opening unit door	Y	
Detachable probe in a bottle filled with a thermal buffer, like bio-safe glycol, which more closely reflects vaccine temperatures.	Y	
Alarm to alert out-of-range temperatures = visual alarm on screen	Y	
Accuracy within $\pm 1.0^{\circ}\text{F}$ OR $\pm 0.5^{\circ}\text{C}$	Y	
Low battery indicator	Y	
Continuous monitoring and recording capabilities to track and record temperatures over time	Y	
Display of current, as well as minimum and maximum temperatures, which indicate the coldest and warmest temperatures recorded for the past 60 days	Y	
Devices must have current and valid ISO 17025 certificate of calibration	Y	

Additional Features
Local Audio Alarm repeats sound every 4 Seconds
CDC Inspections marks, simply pushing a button
Automatically email & FTP transfer report information, including IIS & VTMC (Vaccine Temperature Monitoring Cloud)
Data Logger Profiles, eliminates entry errors. Files can be easily transferred via email
Data Logger configuration – password protection

AIM does not endorse any products featured on the AIM Virtual Exhibit Hall Webinars. The information displayed is provided by vendors and not verified by AIM.

Find more information and event archives at:
<http://www.immunizationmanagers.org/?page=VirtualExhibitHall>

Purchasing Information
Price Range: \$145 for complete monitoring kit
Calibration/ Operating Life: ISO17025:2005 2 year calibration per manufacture, operating life up to 10 years
Bulk ordering options: Yes
Purchasing information: www.vfcdataloggers.com (888) 311-0636
Availability on state contracts: yes
Availability on third party contracts: yes

Technical Support
End User Initial Support Power point & recorded video for setup instructions Webinar training sessions
End User On-going Support (888) 311-0636 Toll FREE TAC 6:30AM to 5:00PM PST
Immunization Program Initial Support Roll out planning, pre-configuration, shipping instructions, administrative training
Immunization Program On-going Support
Existing IT problems: NONE
Change the date/ time/ or C to F Date/time is synced every time logger is docked.. Changing F/C is configurable.

AIM does not endorse any products featured on the AIM Virtual Exhibit Hall Webinars. The information displayed is provided by vendors and not verified by AIM.

Find more information and event archives at:
<http://www.immunizationmanagers.org/?page=VirtualExhibitHall>